



ŞİMŞEK Laborteknik

SAĞ. VE LAB. CİH. SAN. TİC. LTD. ŞTİ.

Merkez Sanayi Sitesi 1345 Sok. No: 22-24 İvedik Org. San. - ANKARA - TÜRKİYE • Tel:+90 312 395 24 84 - 395 24 85
Ulus Vergi Dairesi: 813 006 9397 • Ticaret Sicil No: 114898 • Oda Sicil No: 3/1162 • Faks: 0312 395 24 86
Mersis No: 0813006939700017 • www.simseklabortechnik.com.tr • e-mail: info@simseklabortechnik.com.tr

ANALİZ ELEKLERİ – Laboratuvar Tipi (Paslanmaz)

Şimşek Laborteknik Marka

Elek Ölçüsü Mesh No

32 µm	450 Mesh
38 µm	400 Mesh
45 µm	325 Mesh
53 µm	270 Mesh
63 µm	230 Mesh
75 µm	200 Mesh
90 µm	170 Mesh
106 µm	140 Mesh
125 µm	120 Mesh
150 µm	100 Mesh
180 µm	80 Mesh
212 µm	70 Mesh
250 µm	60 Mesh
300 µm	50 Mesh
355 µm	45 Mesh
425 µm	40 Mesh
500 µm	35 Mesh
600 µm	30 Mesh
710 µm	25 Mesh
850 µm	20 Mesh
1.00 mm	18 Mesh
1.18 mm	16 Mesh
1.40 mm	14 Mesh
1.70 mm	12 Mesh
2.00 mm	10 Mesh
2.36 mm	8 Mesh
2.80 mm	7 Mesh
3.35 mm	6 Mesh
4.00 mm	5 Mesh
4.75 mm	4 Mesh
5,17 mm	1/8"
6.30 mm	1/4"
6.70 mm	0.265"
8.00 mm	5/16"
9.50 mm	3/8"
11.20 mm	7/16"
12.50 mm	1/2"
13.20 mm	0.530"
16.00 mm	5/8"
19.00 mm	3/4"
22.40 mm	7/8"
25.00 mm	1"
31.50 mm	1 1/4"
37.50 mm	1 1/2"
45.00 mm	1 3/4"
50.00 mm	2"
53.00 mm	2.12"
63.00 mm	2 1/2"
75.00 mm	3"
90.00 mm	3 1/2"
100.00 mm	4"



ÜRETİM YAPILAN KASNAK ÖLÇÜLERİ: Çap: 200mm / 300mm / 400mm / 450mm
300x300mm Çubuk Kare Elek // 400x600mm Elek + Toplama Kapları (Kapaklı)